

Product group: 661      Product number: 777932

The most basic test for bacterial contamination of a water supply is the test for total coliform bacteria. Total coliform counts give a general indication of the sanitary condition of a water supply.

For consistent results, an incubator should be used to culture the samples and the correct temperature.

(UV lamp supplied separately)



Product information

Colilert test is EPA approved method for total coliform and E.Coli in drinking water and E.Coli detection in wastewater.

Test Parameters: Coliform and E.Coli bacteria tests

Supplied with box containing pack of 25 Coli test bags and 25 ampules of growth media.

Used with:

- 1. POTABLE WATER TEST KIT (778419)
- 2. BALLAST CHECK 2 TEST KIT (777848)

Features

- Simple presence/absence test
- Easy to carry out on board ship
- Can be to test for Potable water and ballast water
- For Ballast water, dilute the sample by 1 part sample and 9 part deionised water (10 ml sample + 90 ml demin water = 100 ml sample)

Benefits

- Convenient method for measuring and monitoring bacteria and microbial activity in water,
- High quality and consistently reliable results

Specification

General

|                                                |    |
|------------------------------------------------|----|
| Invent Hazard Material (IMO/IU) classification | NA |
|------------------------------------------------|----|

Documents

[SDoC and MD for IHM](#)

Related products

Is accessory to

778419  
POTABLE WATER TEST KIT  
777848  
BALLAST CHECK 2 TEST KIT

Is frequently bought together with

765018  
ENVIROCLEAN 25 LTR  
777934  
CHLORINE DPD NO 1 (FREE) 250 TESTS  
778414  
DILUTION POTS (40 PCS)  
778483  
1 ML STERILE SYRINGES (30)  
778413  
BACTERIA PLATES (40 PCS)

This page is printed from:

<https://www.w.wilhelmsen.com/product-catalogue/products/marine-chemicals/test-kits-and-reagents/ballast-w ater-testing/e.colicoliform-test-set25-tests/>