

# E. COLI/COLIFORM (25 TESTS)

Product group: **661** Product number: **777932**

The most basic test for bacterial contamination of a water supply is the test for total coliform bacteria. Total coliform counts give a general indication of the sanitary condition of a water supply.



For consistent results, an incubator should be used to culture the samples and the correct temperature.

(UV lamp supplied separately)

## Product information

Colilert test is EPA approved method for total coliform and E.Coli in drinking water and E.Coli detection in wastewater.

Test Parameters: Coliform and E.Coli bacteria tests

Supplied with box containing pack of 25 Coli test bags and 25 ampules of growth media.

Used with:

POTABLE WATER TEST KIT (778419)  
BALLAST CHECK 2 TEST KIT (777848)

### Features

- Simple presence/absence test
- Easy to carry out on board ship
- Can be to test for Potable water and ballast water
- For Ballast water, dilute the sample by 1 part sample and 9 part deionised water (10 ml sample + 90 ml demin water = 100 ml sample)

### Benefits

- Convenient method for measuring and monitoring bacteria and microbial activity in water,
- High quality and consistently reliable results

## Specification

### General

Invent Hazard Material (IMO/EU) classification	NA
--	----

## Documents

[SDoC and MD for IHM](#)

## Related products

Is accessory to

**778419**  
POTABLE WATER TEST KIT  
**777848**  
BALLAST CHECK 2 TEST KIT

Is frequently bought together with

**778483**  
1 ML STERILE SYRINGES (30)  
**765018**  
ENVIROCLEAN 25 LTR  
**777934**  
CHLORINE DPD NO 1 (FREE) 250 TESTS  
**778414**  
DILUTION POTS (40 PCS)  
**778413**  
BACTERIA PLATES (40 PCS)

This page is printed from

<https://www.wilhelmsen.com/product-catalogue/products/marine-chemicals/test-kits-and-reagents/ballast-water-testing/e.coliform-test-set25-tests>