

Last updated: 02/11/2023 02:01:07

# DE-IONISED WATER 1 LTR

Product group: **661** Product number: **766414**

Consumables: 1 lit bottle of De-ionised water.

Used to rinse pH electrodes and Conductivity probes.

It can also be used with MAXIDRIECT MD600 as a blank or for mixing of reagents for tests.



## Product information

Deionized (DI) water is water that has been treated to remove all ions – typically, that means all of the dissolved mineral salts. It is suitable for many applications, including autoclaves, laboratory testing - mixing of test reagents etc.

Supplied as 1L bottle

### Features

- Easy and portable for direct use

### Benefits

- Convenient way to have clean, uncontaminated water for cleaning and mixing of test reagents

## Specification

### General

|  |     |
|--|-----|
| Invent Hazard Material (IMO/EU) classification | C-9 |
|--|-----|

## Documents

[SDoC and MD for IHM](#)

## Related products

Is accessory to

**607801**  
CONDUCTIVITY METER PEN MODEL

**777079**  
ECOSCAN CONDUCTIVITY METER MO347

**607800**  
PH METER PEN MODEL

**777080**  
ECOSCAN.PH6 METER & ELECTRODE MO34

Is frequently bought together with

**568683**  
STD CONDUCTIVITY SOLN FOR CALIBRATN 100 ML

**568691**  
COND NEUTRALISING SOLN 65 ML

**555631**  
NITRITE NO.2 TABLETS (2 BOTTLES X250)

**571356**  
ROCOR NB LIQUID 25 LTR

**571752**  
NATURAL HANDCLEANER 4X5 LTR