

Product group: 661 Product number: 766414

Consumables: 1 lit bottle of De-ionised water.
Used to rinse pH electrodes and Conductivity probes.
It can also be used with MAXIDRIECT MD600 as a blank or for mixing of reagents for tests.



Product information

Deionized (DI) water is water that has been treated to remove all ions – typically, that means all of the dissolved mineral salts. It is suitable for many applications, including autoclaves, laboratory testing - mixing of test reagents etc.

Supplied as 1L bottle

Features

- Easy and portable for direct use

Benefits

- Convenient way to have clean, uncontaminated water for cleaning and mixing of test reagents

Specification

General

Invent Hazard Material (IMO/EU) classification	C-9
--	-----

Documents

[SDoC and MD for IHM](#)

Related products

Is accessory to

- 607801
CONDUCTIVITY METER PEN MODEL
- 777079
ECOSCAN CONDUCTIVITY METER MO347
- 607800
PH METER PEN MODEL
- 777080
ECOSCAN.PH6 METER & ELECTRODE MO34

Is frequently bought together with

- 568683
STD CONDUCTIVITY SOLN FOR CALIBRATN 100 ML
- 568691
COND NEUTRALISING SOLN 65 ML
- 555631
NITRITE NO.2 TABLETS (2 BOTTLES X 250)
- 571356
ROCOR NB LIQUID 25 LTR
- 571752
NATURAL HANDCLEANER 4X5 LTR

This page is printed from:

<https://www.wilhelmsen.com/product-catalogue/products/marine-chemicals/test-kits-and-reagents/water-test-kit---spares-and-consumables/de-ionised-water-1-ltr/>