

Product group: **661** Product number: **766414**

Consumables: 1 lit bottle of De-ionised water.

Used to rinse pH electrodes and Conductivity probes.

It can also be used with MAXIDRIECT MD600 as a blank or for mixing of reagents for tests.



## Product information

Deionized (DI) water is water that has been treated to remove all ions - typically, that means all of the dissolved mineral salts. It is suitable for many applications, including autoclaves, laboratory testing - mixing of test reagents etc.

Supplied as 1L bottle

### Features

- Easy and portable for direct use

### Benefits

- Convenient way to have clean, uncontaminated water for cleaning and mixing of test reagents

## Specification

### General

Invent Hazard Material (IMO/EU) classification	C-9
--	-----

## Documents

[SDoC and MD for IHM](#)

## Related products

### Is accessory to

**607801**  
CONDUCTIVITY METER PEN MODEL

**777079**  
ECOSCAN CONDUCTIVITY METER MO347

**607800**  
PH METER PEN MODEL

**777080**  
ECOSCAN.PH6 METER & ELECTRODE MO34

### Is frequently bought together with

**568683**  
STD CONDUCTIVITY SOLN FOR CALIBRATN 100 ML

**568691**  
COND NEUTRALISING SOLN 65 ML

**555631**  
NITRITE NO.2 TABLETS (2 BOTTLES X 250)

**571356**  
ROCOR NB LIQUID 25 LTR

**571752**  
NATURAL HANDCLEANER 4X5 LTR

### This page is printed from:

<https://www.wilhelmsen.com/product-catalogue/products/marine-chemicals/test-kits-and-reagents/water-test-kit---spares-and-consumables/de-ionised-water-1-ltr/>